

OLYMPIC CIRCUMSTRAINT™

*The Safe & Secure Way
to Immobilize Newborns*

**Ideal for Circumcisions,
Transfusions, Minor
Surgery, X-rays, Scalp IVs
& Other Procedures**



*Proven over many years,
the Circumstraint device
has become the infant
immobilization standard
for modern nurseries.*

Versatile & Easy-to-use

Comfortable molded shape elevates hips & presents genitalia for procedure

Perfect for x-rays — minimizes movement, no artifacts

Molded head contour exposes scalp for IVs

Comfortable, adjustable straps securely hold infant at elbows & knees

Durable plastic construction ensures years of use

Lightweight & compact design enables easy handling, storing & cleaning



OLYMPIC CIRCUMSTRAINT™

First introduced in 1959, the Olympic Circumstraint device has been purchased by over 11,000 hospitals throughout the world. The compact Circumstraint unit, always a fast, easy means of immobilizing newborns for circumcision, fits inside a standard incubator or bassinet and is useful for a variety of procedures. The Circumstraint unit may be used for minor surgery, blood exchanges, catheterization, transfusions or infusions (including scalp IV's), or whenever immobilization is required for a neonatal procedure. It also works very well for x-rays, as it is radiolucent and holds the infant very securely with minimal movement.

In less than 30 seconds, a nurse can immobilize an infant securely in the correct position with the Circumstraint device. Soft, wide Velcro® straps encircle an infant's elbows and knees, ensuring safe & secure immobilization. The Circumstraint unit's comfortable contoured shape positions the infant, hips elevated, effectively presenting the genitalia for the circumcision procedure. The platform between the infant's legs provides support for a circumcision clamp. Without extraneous pins, clamps or towels, the device snugly and securely immobilizes the infant with the entire torso visible.

The Circumstraint device fits infants from 1800 grams to 4500 grams. Comfortable, adjustable straps adapt to infant length, while Velcro® closures easily adjust to extremity circumference. The device maintains the high standards of nursery asepsis with minimum effort. The Circumstraint unit is constructed of dense, smooth, non-absorbent plastic and wipes clean with liquid disinfectant. Its nylon Velcro® straps dry quickly, and may be laundered, steam or gas autoclaved, or pasteurized. The device is very cost-effective, as it permits a nurse to quickly immobilize an infant in the proper position, speeding procedures and saving valuable time. Exceptionally durable, Circumstraint provides years of service, even with long-term use. Proven over many years, Circumstraint has become the infant immobilization standard for modern nurseries.



One-Year Guarantee: The Olympic Circumstraint device is fully guaranteed for one year.

Specifications & Ordering Information

OLYMPIC CIRCUMSTRAINT™

Lightweight, compact & easy to hang or store in the OR or nursery. Several Circumstraint units can stack nested on a shelf.

Dimensions: 25¼" x 13" x 3½" (64.5 cm x 33 cm x 8.9 cm)

Weight: Less than 2 lbs (0.9kg)

Circumstraint Newborn Immobilizer

Complete with straps.

CAT. NO. 50100 \$265.00

Regular Replacement Strap Set

Set of 8 straps: Four 8" (20.3 cm) arm straps and four 10" (25.4 cm) leg straps.

CAT. NO. 50108 \$49.00

Longer Leg Replacement Strap Set

Set of 8 straps: Four 8" (20.3 cm) arm straps and four 13" (33cm) leg straps.

CAT. NO. 50109 \$49.00

natus
pediatrics

Natus Medical Incorporated
1501 Industrial Road • San Carlos, CA 94070 USA
1-800-303-0306 • 1-650-802-0400
www.natus.com

For information or ordering
Call 1-800-303-0306
(toll-free in US/Canada)