

ATTENTION

USE THESE HEEL LANCETS

**FOR LAB DRAWS.
NOT FOR GLUCOSE!**

**FOR GLUCOSE.
FOR A FEW DROPS OF BLOOD.**



Full-Term

For Babies Weighing
MORE THAN
1,800-2,000 grams
(4-5 lbs)
Mango Orange

Length 2.5mm-3.0mm
Depth 1.0mm-1.25mm



Premie

For Babies Weighing
LESS THAN
1,800-2,000 grams
(4-5 lbs)
Lime Green

Length 1.75mm-2.25mm
Depth .85mm-1.0mm



NeatNick™



SugarPlum™



Or Use Your Hospital's Glucose Lancet.

Hawaii Medical NeatNick

Babies Weighing LESS THAN 1,800 – 2,000 Grams

1030083 Premie Size Lime Green 100 per box

Babies Weighing MORE than 1,800 – 2,000 Grams

1030085 Full-Term Size Mango Orange 100 per box

Hawaii Medical SugarPlum

Micropreemies to Full-Term Babies

1040043 Un-bagged, 1.0mm depth (yellow) 100 per box

1040044 Individually Bagged, 1.0mm depth (yellow) 100 per box

Full-Term Babies and Adults

1040041 Un-bagged, 1.8mm depth (gray) 100 per box

1040042 Individually Bagged, 1.8mm depth (gray) 100 per box



natus.
pediatrics



Natus Medical Incorporated

1501 Industrial Road • San Carlos, CA 94070 USA | 1.800.255.3901 (US/CA) • Worldwide: 1.650.802.0400

www.natus.com

P/N 003401A