

The new NeuroMax CE system expands the XLTEK legacy of clinical utility, durability & portability.

The ultra compact NeuroMax CE system is available in 2 & 4 channel models facilitating complex Neurodiagnostic testing such as:

- Multi-channel EMG Recording
- QEMG
- Multi-channel Nerve Conductions
- Autonomic Nervous System Testing
- Robust Evoked Potential Suite

A frequent comment from owners of all generations of the NeuroMax system: *"It is the most dependable EMG I have ever used"*. This quality is not by accident...*it's by design.*

The NeuroMax CE system:

- Has **no internal moving parts**; no disk drives or other moving components that can wear out or fail over time. It is the only **100% solid-state** EMG system available!
- Is **strong, lightweight** and designed to withstand the test of time. Only the highest grade of impact-resistant PC+ABS resin is used to house NeuroMax CE and its preamplifier.
- Is now **almost 3 lbs (1.5 kg) lighter** than previous generations.
- Features **buttons/keys of extremely high quality**, designed for **20 million operations**.
- Features a keyboard strengthened with metal posts for **improved rigidity & feel**.

For the clinician on-the-go, portability is a top priority. The NeuroMax CE system is the *only portable EMG ready to be used right out of the box.*

The NeuroMax CE system:

- Is **completely self-contained**; no need for cumbersome setup & connection to an external notebook PC.
- Is **ready to use in about 30 seconds**; a typical notebook PC takes at least 2 minutes of start up time before it is ready for use.
- Weighs only 11 lbs (5 kg), making it the **lightest portable EMG system available!**

Contact your local Natus District Sales Manager for more information: 1-800-255-3901

