



Give your babies the comfort of oral care with the gentle foam tip swab applicator

- Foam tip applicator is soft, safe and absorbent for your oral care applicator needs
- Available in two convenient sizes: extra-small 6mm wide and small 9 mm wide
- Intended use is general cleaning and removing bacteria with solution
- Cotton under the foam provides extra absorbency

Babies Enjoy it

Caregivers report that babies enjoy the gentle swabbing action. Oral care takes only a few minutes to administer at touch time, once per shift, or as needed. (Be sure to follow your hospital protocol or request our sample protocol.) In addition, it provides an excellent opportunity for proactive family involvement when family, staff and baby are ready. Parents often are able to perform the oral care protocol and feel involved in their baby's care.

Ordering Information



ITEM #	DESCRIPTION	QUANTITY
004553	Extra-Small 6mm wide, 20mm long and 155mm handle	500 per bag
004554	Small 9mm wide, 20mm long and 155mm handle	500 per bag

PRODUCT SPECIFICATIONS:
Custom medical grade foam over cotton, non-sterile.

Hawaii Medical
OraSwab™

