



Natus Neurology IT Networking Solutions for EMG

Connectivity and workflow solutions designed for you

Natus HL7 Solutions

Natus Neurology works cooperatively with clients to facilitate connectivity solutions to improve the patient experience. The dedicated service team works closely with information technology staff to integrate Natus HL7 compatible systems. Whether a private practice or hospital setting, Natus equipment allows an optimized workflow from acquisition to final report to the hospital Electronic Medical Records System.

Natus Neurology equipment has been successfully integrated into Epic, Meditech, Cerner, McKesson and other Hospital Electronic Medical Records Systems.

Remote Access

Natus EMG systems are compatible with and can be leveraged with the Citrix® XenApp or VMware® vSphere solutions. This allows users to remotely access and review studies, which provides the opportunity to improve cost efficiencies as well as quality and delivery of patient care.

Natus provides the most powerful infrastructure to facilitate collaboration with offsite colleagues through remote access and review.



HIPAA Compliance

Natus EMG systems are compatible with Microsoft Active Directory and Microsoft Bitlocker encryption. Additionally, Natus EMG systems not only allow user-defined, user-based security/configurations for restricted user access to Patient Health Information (PHI), but also provides automatic user activity logging for auditing.

Security can be integrated with Microsoft® Active Directory as well as with Microsoft Bitlocker encryption.

Basic EMR Transfer

Natus EMG systems allow for a fast, easy report transfer option with the Basic EMR Transfer feature in its software. It allows finalized reports to be saved to a centralized folder for EMR transfer via basic scripts. After simple setup of the folder location, just click the transfer icon and confirm to efficiently send pdf reports, streamlining your workflow.

Simplify the process of transferring finalized pdf reports to your EMR with the Basic EMR Transfer feature.



Service

Excellence in customer and technical service is a core value of Natus Neurology.

Essential Support Elements:

- Accessible and effective Technical Support
- Definitive technical documentation and knowledgeable installation teams
- Excellent replacement unit or spare part availability
- Extended warranty and service coverage programs
- Comprehensive, flexible customer training courses

Supplies

Natus Neurology offers a full range of neurodiagnostic supplies, featuring:

- Dedicated and knowledgeable customer service
- Streamlined order processing to save you time
- Convenient online portal (US Customers only)
[Natus Neuro Store – www.natusneurostore.com](http://www.natusneurostore.com)

To learn more about Natus Neurology Service Programs or our full line of Supplies, contact your local distributor or sales representative.

US Customers Call: **1-800-356-0007**

International Customers Call: **+1-608-829-8500**



Natus Neurology Incorporated

3150 Pleasant View Road
Middleton, WI 53562 USA
Tel: 1-800-356-0007
1-608-829-8500
Fax: 1-608-829-8709
www.natus.com

© 2015 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates. Citrix is a registered trademark of Citrix Systems, Inc. VMware is a registered trademark of VMware, Inc.