

# OLYMPIC CFM™ 6000

## Infant aEEG Cerebral Function Monitor

# ACCESSORIES

### LOW IMPEDANCE NEEDLE ELECTRODES

62056 / 62057



Protective  
Cap Feature

### NEONATAL HYDROGEL SENSORS

62050 / 62051



### NUPREP GEL

102566N



### THERMAL PAPER

62066



### ORDERING INFORMATION

Part #	Item Description
62056	Disposable Low Impedance Needle Electrodes - 12mm-29G. 36 inch long lead wires terminating in female 1.5mm touch proof connectors. Lead wire colors: red, yellow, and black. 3 needles/pouch (box of 30 needles).
62057	Disposable Low Impedance Needle Electrodes - 12mm-29G. 36 inch long lead wires terminating in female 1.5mm touch proof connectors. Lead wire colors: red, yellow, and black. 3 needles/pouch (box of 150 needles).
62050	Disposable Neonatal Hydrogel Sensors - 24 inch lead wires terminating in female 1.5mm touch proof connectors. Lead wire colors: red, yellow, and black. 3 sensors/pouch (box of 60 sensors).
62051	Disposable Neonatal Hydrogel Sensors - 24 inch lead wires terminating in female 1.5mm touch proof connectors. Lead wire colors: red, yellow, and black. 3 sensors/pouch (box of 300 sensors).
102566N	NuPrep Skin Preparation Gel - 4oz Tubes (box of 3).
62066	Thermal Paper - 50mm wide roll.

**natus**  
pediatrics

Natus Medical Incorporated  
1501 Industrial Road • San Carlos, CA 94070 USA  
1-800-303-0306 (US/CA) • Worldwide: +1-650-802-0400  
www.natus.com

DOC-005217A