

BRAINZ

ACCESSORIES

NEONATAL NEEDLE ELECTRODES

004833*



*Shown with CZA00012



Protective
Cap Feature

NEONATAL HYDROGEL SENSORS

CZA00037*



*Shown with CZA00012



New
Zip-Lock Pouch!

SKIN MARKERS

CZN00011



WRAP HATS

CZN00012



POSITIONING STRIPS

OBM00047



NUPREP GEL

102566N



ORDERING INFORMATION

Part # **Item Description**

004833	Disposable Subdermal Needle Electrodes - 12mm-29G. 18 inch long lead wires terminating in female 1.5mm touch proof connectors. Lead wire colors: 2 x purple and 2 x black. 4 needles/pouch (box of 24 needles).
CZA00037	Disposable Neonatal Hydrogel Sensors - 12 inch lead wires terminating in female 1.5mm touch proof connectors. Lead wire colors: 5 x white. 5 sensors/pouch (box of 60 sensors) in re-sealable zip-lock pouch.
CZA00012	Sensor Adapter Set with positioning guide (box of 10)
CZN00012	Wrap Hats (bag of 10)
CZN00011	Skin Markers (box of 10)
102566N	NuPrep Skin Preparation Gel - 4oz Tubes (box of 3)
OBM00047	Positioning Strips - Term and Pre-Term, pack of 20 (10 of each)

natus
pediatrics

Natus Medical Incorporated
1501 Industrial Road • San Carlos, CA 94070 USA
1-800-303-0306 (US/CA) • Worldwide: +1-650-802-0400
www.natus.com

DOC-005218B